



## **Friesian Horse Society, Inc. (FHS)** **Instructions for Microchip Implantation**

---

Please provide the following instructions to your veterinarian:

The microchip is implanted in the nuchal ligament on the left side of the neck perpendicular to the point of the shoulder.  
After implantation the chip should be scanned and verified. Please make sure to record the number and place the sticker on the Microchip form.

*If there are any questions please contact the FHS office!*

**FHS - Friesian Horse Society Inc.**  
**17670 Pioneer Trail**  
**Plattsmouth, NE 68048**  
E-mail: [info@friesianhorsesociety.com](mailto:info@friesianhorsesociety.com)

**Phone: 402-298-7347**  
**FHS Fax: 402-298-7569**